

**NEONATAL AND  
PAEDIATRIC RANGE  
SPECIFICALLY DESIGNED  
FOR THEIR ANATOMY**



# AVANOS\* NEONATAL AND PAEDIA

AT AVANOS WE UNDERSTAND THE NEONATAL AND PAEDIATRIC NEEDS. OUR SOLUTIONS BRING YOUR PATIENT SAFETY AND COMFORT WITH A FIRM COMMITMENT TO INFECTION CONTROL

## AVANOS\* MICROCUFF\* PAEDIATRIC ENDOTRACHEAL TUBES


A NEW STANDARD FOR PAEDIATRIC AIRWAY MANAGEMENT

### THE PERFECT FIT FOR PAEDIATRIC ANATOMY


- Cuff placement in the trachea and not in the pressure-sensitive larynx
- Intubation depth markings to ensure optimal tube placement
- Microthin cuff material to seal at lower pressures<sup>1</sup> than conventional tubes reducing the risk to mucosal tissue<sup>2</sup>
- Withstand burst pressures of >800 cm H<sub>2</sub>O<sup>3</sup>

"The intubation depth marks of the MICROCUFF\* paediatric tracheal tube allowed the safe placement of a cuffed tracheal tube in children from a wide age range" - Weiss, et al. Br J Anaesthesia 2005


### MICROCUFF\* tubes seal at a lower pressure than conventional paediatric tubes<sup>1</sup>



Median cuff pressure to seal the trachea in children aged 2-4 (n=4x20 patients, ID 4.0mm). Sealing pressure assessed by auscultation within 5 minutes after intubation.



In an unstandardized market of unequal paediatric tubes, the MICROCUFF\* tube is designed for paediatric airway




# ATRIC SOLUTIONS

Low dead-space adapter  
on sizes  
3.0 - 4.5mm



## CONFIDENCE IN A SEALED AIRWAY

- Low rate of re-intubation<sup>4</sup>
  - Reduces need to replace oversized tracheal tubes
  - Less patient trauma, time and material costs
- Sealed airway allows for use of minimal and low flow anesthesia
- Sealing with a cuff compensates for different sized and shaped airways



## IMPROVED PATIENT CARE AND COMFORT<sup>1</sup>

- Positive pressure ventilation with a sealed airway, providing constant and efficient minute ventilation
- Sealed airway ensures reliable end-tidal CO<sub>2</sub> lung function and oxygen consumption monitoring
- Reduced risk of aspiration of blood and secretions

## AVANOS\* CLOSED SUCTION SYSTEMS FOR NEONATES/PAEDIATRICS

GOLD STANDARD IN CLOSED SUCTION SYSTEMS

Closed suctioning has been extensively studied and has evident advantages over open suctioning for both patients and their caregivers.

- Smaller degree of oxygen saturation fall and fewer incidences of desaturation<sup>5</sup>
- Reduced incidences of bradycardia and hypoxia associated with suctioning<sup>6</sup>
- Standardize the level of care within a ward
- Safely used by one person, reducing the need for two or more caregivers often required for open suctioning procedures<sup>7</sup>
- Perceivably easier, less time-consuming, and better tolerated<sup>7</sup> by small premature infants requiring mechanical ventilation for one or more weeks



**Compared to open suctioning, no increase in the rate of<sup>7</sup>:**

- bacterial airway colonization,
- frequency of endotracheal suctioning and reintubation,
- duration of mechanical ventilation,
- length of hospitalization,
- incidence of nosocomial pneumonia, nosocomial BSI, severity of BPD, and neonatal mortality.<sup>7</sup>

**The unique AVANOS\* Multi-Access Catheter is specially designed to**

- Access the neonatal airway
- Deliver surfactant without disconnection from the ventilator leading to better distribution of surfactant and almost no wastage


**Key features:**

- Numerous universal colour stripes
- Sizes 5 to 12 French  
More flexibility for matching the appropriate sized suction catheter to a size 3.5 endotracheal tube
- Noticeably softer catheter
- DEHP-free - All neonatal and paediatric catheters are DEHP-free



## **AVANOS\* ORAL CARE SOLUTIONS FOR NEONATES/PAEDIATRICS**

Avanos also offers a Neonatal and Paediatric Oral Care range of small, soft, atraumatic catheters, designed to clean the oral cavity effectively in smaller patients whilst keeping the system closed.




## AVANOS\* MICROCUFF\* PAEDIATRIC ENDOTRACHEAL TUBES

### AVANOS\* MICROCUFF\*

Endotracheal Tube, Paediatric Oral/Nasal Magill

REFERENCE	TUBE SIZE (INNER DIAMETER)
<b>35111</b>	3.0 mm
<b>35112</b>	3.5 mm
<b>35113</b>	4.0 mm
<b>35114</b>	4.5 mm
<b>35115</b>	5.0 mm
<b>35116</b>	5.5 mm
<b>35117</b>	6.0 mm
<b>35118</b>	6.5 mm
<b>35119</b>	7.0 mm

Packaging of each item: 10/dispenser



### AVANOS\* MICROCUFF\*

Endotracheal Tube, Paediatric Oral Curved

REFERENCE	TUBE SIZE (INNER DIAMETER)
<b>35161</b>	3.0 mm
<b>35162</b>	3.5 mm
<b>35163</b>	4.0 mm
<b>35164</b>	4.5 mm
<b>35165</b>	5.0 mm
<b>35166</b>	5.5 mm
<b>35167</b>	6.0 mm
<b>35168</b>	6.5 mm
<b>35169</b>	7.0 mm

Packaging of each item: 10/dispenser



## AVANOS\* CLOSED SUCTION SYSTEMS FOR NEONATES/PAEDIATRICS

### AVANOS\* Closed Suction Systems

for Neonates/Paediatrics - "Y"

REFERENCE	OUTER DIAMETER	LENGTH
<b>195-5</b>	5 Fr/1.6mm	30.5cm
<b>196-5</b>	6 Fr/2mm	30.5cm
<b>197-5</b>	7 Fr/2.3mm	30.5cm
<b>198-5</b>	8 Fr/2.6mm	30.5cm
<b>1910-5</b>	10 Fr/3.3mm	40.5cm
<b>1912-5</b>	12 Fr/4mm	40.5cm

Packaging of each item: 20/case



### AVANOS\* Closed Suction Systems

for Neonates/Paediatrics - Elbows

REFERENCE	OUTER DIAMETER	LENGTH
<b>206-5</b>	6 Fr/2mm	30.5cm
<b>207-5</b>	7 Fr/2.3mm	35.5cm
<b>208-5</b>	8 Fr/2.6mm	35.5cm
<b>210-5</b>	10 Fr/3.3mm	40.5cm

Packaging of each item: 20/case



### Multi-Access Catheter

REFERENCE	OUTER DIAMETER	LENGTH	DESCRIPTION
<b>1900</b>	5 Fr/1.6mm	30.5cm	MAC catheter adapter: 2.5 mm, 3 mm, 3.5 mm
<b>1920</b>	5 Fr/1.6mm	30.5cm	MAC catheter adapter: 2 mm
<b>1940</b>	5 Fr/1.6mm	30.5cm	MAC catheter adapter: 4 mm

Packaging of each item: 5/case



## AVANOS\* ORAL CARE SOLUTIONS FOR NEONATES/PAEDIATRICS

REFERENCE	DESCRIPTION
1221	5 oral care suction catheters: 6 Fr - 2 mm
1222	5 oral care suction catheters: 8 Fr - 2.6 mm
1223	5 oral care suction catheters: 10 Fr - 3.3 mm

Packaging of each item: 250/case

All codes in **bold** indicate a **sterile** product.

**References** 1. Dullenkopf A, Schmitz A, Gerber A, Weiss M. Tracheal sealing characteristics of paediatric cuffed tracheal tubes. *Paediatric Anesthesia*. 2004; 14:825-830. 2. Seegobin RD, van Hasselt GL. Endotracheal cuff pressure and tracheal mucosal blood flow: endoscopic study of effects of four large volume cuffs. *British Medical Journal*. 1984 March;228:965-968. 3. Data on file. Roswell, GA, KCWW. 4. Dullenkopf A, Gerber AC, Weiss M. Fit and seal characteristics of a new paediatric tracheal tube with a high volume-low pressure polyurethane cuff. *Acta Anaesthesiol Scand*. 2005;49:232-237. 5. Tan, A.M., Gomez, J.M., Matthews, J., Williams, M., Paratz, J., Rajadurai VS. (2005). Closed versus partially ventilated endotracheal suction in extremely preterm neonates: physiologic consequences. 6. Wright, J. & Fraser Askin, D. (1996). Closed-suctioning procedure in neonates. 7. Cordero, L., Sananes, M. and Avers, L.W. (2000). Comparison of a closed (Trach Care MAC) with an open endotracheal suction system in small premature infants.

---

**For more information, please contact  
your Avanos Customer Service Team:**

**From UK & IE:**  
**tel: 0800 917 65 85 – fax: 0800 169 02 35**  
**customerservice.uk.ie@avanos.com**

**From other countries:**  
**tel: +32 2 700 68 51 – fax: +32 2 711 26 91**  
**customerservice.export@avanos.com**

**www.avanos.co.uk**

**AVANOS**