

**PHILIPS**

RESPIRONICS

**Wisp**

Pediatric nasal mask



**Rediscover  
dreams**



## Meet Jacky Giraffe and friends

Jacky Giraffe is a fun character based on the look of the mask itself. Jacky stars in an animated cartoon that children and their parents can watch on DreamMapper or the Philips Respironics YouTube channel. Kids will see Jacky exploring a fantasy world with Sami the Seal and Tucker the Turtle—promoting acceptance of the mask.



## Fitting guide for Wisp pediatric patients

Achieving the right fit with Wisp pediatric is easy. Simply follow the instructions below.



1. The mask cushion should fit the width of the nose without blocking the nostrils.
2. Pulling the headgear over the child's head, gently hold the cushion over the nose.
3. Reconnect the clips by pushing them into the mask frame.
4. Connect the Wisp pediatric tubing and elbow swivel to the front of the mask.

*Three cushion sizes (SCS, SCM, SCL) are available. Please read the Instructions for Use for additional fitting information.*



## Minimal contact, **maximum comfort**

Get the comfort and performance of the Wisp minimal contact nasal mask—designed specifically for children.

Over the years, adult masks have become smaller, lighter, and less intrusive, while pediatric masks have largely remained the bulky masks of the past. The Wisp pediatric nasal mask is truly designed for children. Wisp pediatric applies a child-friendly fabric pattern, a modified cushion curvature, and family support tools to provide a positive experience for your smallest patients.

Now, with the availability of Wisp pediatric, Wisp youth, and Wisp minimal contact nasal masks, the Wisp platform serves infants to adults, providing continuity and consistency in their therapy experience.

# Unique design

for unique patient needs

**Wisp pediatric evolves the proven Wisp nasal mask design for the pediatric population.**


- The cushion, frame, and headgear are designed to fit the specific sizing and bone structure of pediatric patients.
- The three cushion sizes are designed specifically for young children.
- A Leak Correction Dial lets parents adjust the mask without removing it or waking the child.
- An open line of sight allows children to comfortably wear the mask during their normal daily activities.



**Optional  
tube wrap with tubing  
management  
loop**

Enhances the giraffe print theme while adding functionality

The tube wrap with tubing management loop enables the tube to be positioned over the head and away from the front of the body.



**360° rotating elbow  
for freedom of movement**  
Features a quick release  
design

**Cushion design,  
featuring auto seal groove  
technology**  
Delivers minimal contact  
and a snug seal

**Click-fit headgear  
clips with adjustable straps**  
Provides for quick and easy  
fitting and adjustment

**Leak Correction Dial**  
Enables less intrusive  
mask adjustment

**The innovative Leak Correction Dial on the Wisp pediatric mask makes fixing small leaks easy.**

Turn it to the right or left if leaks occur while the child is sleeping to readjust the fit of the cushion on the face. This widens the angle of the front of the frame, gently nudging the cushion closer to the face.

# Simplified ordering

for increased efficiency

**Wisp pediatric features a streamlined product offering to help make it easier to stock and supply to your customers.**

Give your entire pediatric population the comfort they need with Wisp pediatric and Wisp youth.

## Wisp pediatric nasal mask

	Part number
Fit pack with 3 cushion sizes	1104953

## Replacement parts

Headgear	1104973
Elbow/tube with cover	1104977
Frame	1104972
Small cushions (SCS)	1104969
Medium cushions (SCM)	1104970
Large cushions (SCL)	1104971

## Wisp youth nasal mask

	Part number
Fit pack with 3 cushion sizes	1109298

For more information, call

**1-800-345-6443**



© 2016 Koninklijke Philips N.V. All rights reserved.  
Specifications are subject to change without notice.

[www.philips.com/respironics](http://www.philips.com/respironics)

CAUTION: US federal law restricts these devices  
to sale by or on the order of a physician.

RRDPGH EB 4/1/16 MCI 4107201 PN 1127222  
1010 Murry Ridge Lane, Murrysville, PA 15668  
800-345-6443 • 724-387-4000