

## GO-PAP™

### Emergency Disposable CPAP with Integrated Nebulization

- Runs at 10 LPM of oxygen, and generates up to 90 LPM of flow.
- Has an approximate 30% FiO<sub>2</sub>.
- Includes 3 options for PEEP of 5, 7.5, and 10 cm H<sub>2</sub>O.
- Completely disposable system.
- Run time of approximately 40 minutes.
- Integrated nebulization without interrupting the delivered flow to the patient.
- Nebulize off of one tank utilizing the Neb-Connect™.

## Neb-Connect™



CPAP & NEBULIZE  
Using **ONE**  
Tank!

### GO-PAP™ System with BiTrac™ ED Mask and Neb-Connect™- nebulizer canister included

Part No.	Headgear	Size	Case/Qty.
313-4604NC	CAP	Adult Medium	10
313-4605NC	CAP	Adult Large	10
313-4602NC	Standard	Adult Medium	10
313-4603NC	Standard	Adult Large	10

### Neb-Connect™

Part No.	Recommended Flow Rate	Case/Qty.
313-011111	6-8 LPM	30

### GO-PAP™ System with BiTrac™ ED Mask - nebulizer canister included

Part No.	Headgear	Size	Case/Qty.
313-4601N	Standard	Adult Small	10
313-4602N	Standard	Adult Medium	10
313-4603N	Standard	Adult Large	10

### GO-PAP™ System with BiTrac™ ED Mask - no nebulizer canister included

Part No.	Headgear	Size	Case/Qty.
313-4601	Standard	Adult Small	10
313-4602	Standard	Adult Medium	10
313-4603	Standard	Adult Large	10



CAP Headgear



Standard Headgear