

# Hamilton accessories

## CONTENT

- 2 | Breathing circuit set
- 3 | Breathing circuit set for HAMILTON-H900
- 3 | Flowsensor
- 4 | Expiratory valve & membrane
- 5 | CO<sub>2</sub>
- 6 | SpO<sub>2</sub>
- 6 | Nebulising
- 7 | Masks
- 8 | Transpulmonary pressure
- 8 | Intellicuff
- 8 | Protective sleeve
- 8 | nCPAP pressure line





260087 / 260094 / 260145 / 260144



260207 / 260240



260128 / 260127 / 260167 / 260168



260257 / 260257



260180 / 260182



260170 / 260169




## Breathing circuit set

|        |                                                                                                                                                                                              |   |           | C1/T1 | C3/C6 | G5/S1 | MR1 |
|--------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---|-----------|-------|-------|-------|-----|
| 260087 | <b>Coaxial + flow 180 cm.</b> Pre-assembled set with elbow adapter and flow sensor (without expiratory valve)                                                                                | ⊗ | Box of 20 | X     |       |       | X   |
| 260094 | <b>Coaxial + flow 240 cm.</b> Pre-assembled set with elbow adapter and flow sensor (without expiratory valve)                                                                                | ⊗ | Box of 10 | X     |       |       | X   |
| 260145 | <b>Coaxial + flow 300 cm.</b> Pre-assembled set with elbow adapter and flow sensor (without expiratory valve)                                                                                | ⊗ | Box of 10 | X     |       |       | X   |
| 260144 | <b>Coaxial + flow 480 cm.</b> Pre-assembled set with elbow adapter and flow sensor (without expiratory valve)                                                                                | ⊗ | Box of 8  | X     |       |       | X   |
| 260207 | <b>Coaxial + flow 180 cm.</b> Pre-assembled set with expandable expiratory limb and flow sensor                                                                                              | ⊗ | Box of 20 | X     | X     | X     | X   |
| 260240 | <b>Coaxial + flow 240 cm.</b> Pre-assembled set with expandable expiratory limb and flow sensor                                                                                              | ⊗ | Box of 10 | X     | X     | X     | X   |
| 260128 | <b>Coaxial + flow + expir. valve 180 cm.</b> Pre-assembled set with elbow adapter, flow sensor and expiratory valve                                                                          | ⊗ | Box of 20 | X     |       |       | X   |
| 260184 | <b>Coaxial + flow + expir. valve 180 cm.</b> Pre-assembled set with expandable expiratory limb, flow sensor and expiratory valve                                                             | ⊗ | Box of 20 |       |       | X     |     |
| 260127 | <b>Coaxial + flow + expir. valve 240 cm.</b> Pre-assembled set with elbow adapter, flow sensor and expiratory valve                                                                          | ⊗ | Box of 10 | X     |       |       | X   |
| 260167 | <b>Coaxial + flow + expir. valve 300 cm.</b> Pre-assembled set with elbow adapter, flow sensor and expiratory valve                                                                          | ⊗ | Box of 10 | X     |       |       | X   |
| 260168 | <b>Coaxial + flow + expir. valve 480 cm.</b> Pre-assembled set with elbow adapter, flow sensor and expiratory valve                                                                          | ⊗ | Box of 8  | X     |       |       | X   |
| 260257 | <b>Coaxial + flow + expir. valve 180 cm.</b> Pre-assembled set with expandable expiratory limb, flow sensor and expiratory valve                                                             | ⊗ | Box of 20 |       | X     |       |     |
| 260180 | <b>Neonatal + dual limb 150 cm.</b> Pre-assembled set with Y-piece, flow sensor, flow sensor calibration adapter, pressure line with T-piece connectors                                      | ⊗ | Box of 20 | X     |       |       | X   |
| 260182 | <b>Neonatal + dual limb 300 cm.</b> Pre-assembled set with Y-piece, flow sensor, flow sensor calibration adapter, pressure line with T-piece connectors                                      | ⊗ | Box of 10 | X     |       |       | X   |
| 260170 | <b>Neonatal + dual limb + expir. valve 150 cm.</b> Pre-assembled set with Y-piece, flow sensor, flow sensor calibration adapter, pressure line with T-piece connectors                       | ⊗ | Box of 20 | X     |       |       | X   |
| 260169 | <b>Neonatal + dual limb + expir. valve 300 cm.</b> Pre-assembled set with swivel-Y-piece, flow sensor, flow sensor calibration adapter, pressure line with T-piece connectors with luer lock | ⊗ | Box of 10 | X     |       |       | X   |


⊗ Single Use

↻ Reusable


AC Autoclavable




260161



260186



260185



260187



281637 / 282049



155500 / 260179



279937



279964

## Breathing circuit set for HAMILTON-H900

|        |                                                                                                                                                                                                       |   |           | C1/T1 | C3/C6 | G5/S1 | MR1 |
|--------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---|-----------|-------|-------|-------|-----|
| 260161 | <b>Adult heated dual limb breathing set.</b> Pre-assembled, heated breathing circuit with water chamber, temperature probe and all-in-one connector. For patients over 10 kg                          | ⊗ | Box of 15 | X     | X     | X     |     |
| 260186 | <b>Adult heated single limb breathing set.</b> Pre-assembled, heated breathing circuit with water chamber, temperature probe and all-in-one connector. For patients over 10 kg                        | ⊗ | Box of 15 | X     | X     | X     |     |
| 260185 | <b>Infant heated dual limb breathing set.</b> Pre-assembled, heated breathing circuit with water chamber, incubator extension, temperature probe and all-in-one connector. For patients up to 15 kg   | ⊗ | Box of 15 | X     | X     | X     |     |
| 260187 | <b>Infant heated single limb breathing set.</b> Pre-assembled, heated breathing circuit with water chamber, incubator extension, temperature probe and all-in-one connector. For patients up to 15 kg | ⊗ | Box of 15 | X     | X     | X     |     |

## Flowsensor

|        |                                                                   |   |           |   |   |   |   |
|--------|-------------------------------------------------------------------|---|-----------|---|---|---|---|
| 281637 | <b>Flowsensor. Adult/Ped. 188 cm.</b> 15OD x 15ID/22OD connectors | ⊗ | Box of 10 | X | X | X | X |
| 282049 | <b>Flowsensor. Adult/Ped. 260 cm.</b> 15M x 15F/22M connectors    | ⊗ | Box of 10 | X | X | X | X |
| 155500 | <b>Flowsensor. Infant/Neo. 188 cm.</b> 15M x 15F (patient side)   | ⊗ | Box of 10 | X | X | X | X |
| 260179 | <b>Flowsensor. Infant/Neo. 300 cm.</b> 15M x 15F (patient side)   | ⊗ | Box of 10 | X | X | X | X |
| 279937 | <b>Calibration adapter for flow sensor. Adult/Ped.</b>            | ⊗ | Box of 10 | X | X | X | X |
| 279964 | <b>Calibration adapter for flow sensor. Infant/Neo.</b>           | ⊗ | Box of 10 | X | X | X | X |





## Expiratory valve & membrane

|           |                                                                    |   |           | C1/T1 | C3/C6 | G5/S1 | MR1 |
|-----------|--------------------------------------------------------------------|---|-----------|-------|-------|-------|-----|
| 950158    | <b>Expiratory valve set.</b> With cover and membrane. Single use   | ⊗ | Box of 10 |       |       | X     |     |
| 160176    | <b>Expiratory valve set.</b> With housing and membrane. Single use | ⊗ | Box of 10 |       | X     |       |     |
| 161186    | <b>Expiratory valve set.</b> With cover and membrane. Single use   | ⊗ | Box of 10 | X     |       |       | X   |
| 161189    | <b>Expiratory valve set.</b> With cover and membrane. Single use   | ⊗ | Box of 10 | X     |       |       | X   |
| 151233    | <b>Expiratory valve membrane.</b> For use with 151228              | ⊗ | Box of 5  |       |       | X     |     |
| 160245    | <b>Expiratory valve set.</b> With cover and membrane               | ⊗ | 1 pc      |       | X     |       |     |
| 160500    | <b>Expiratory valve membrane.</b> For use with 160245              | ⊗ | Box of 5  |       | X     |       |     |
| 161175    | <b>Expiratory valve set.</b> With cover and membrane.              | ⊗ | 1 pc      | X     |       |       | X   |
| 161188    | <b>Expiratory valve set.</b> With cover and membrane               | ⊗ | 1 pc      | X     |       |       | X   |
| 161390    | <b>Expiratory valve membrane.</b> For use with 161175 + 161188     | ⊗ | Box of 5  | X     |       |       | X   |
| 151972    | <b>Expiratory valve set.</b> With cover and membrane               | ⊗ | 1 pc      |       |       | X     |     |
| 100074247 | <b>Adapter</b> for Expiratory Valve                                | ⊗ | Box of 30 | X     | X     |       |     |



Single Use



Reusable



Autoclavable



281928



8000.0363



8000.0364



282157



281718



8000.0260.01



8000.0261.01



8000.0262.01



281803

## CO2

|              |                                                                                                                                |  |           | C1/T1 | C3/C6 | G5/S1 | MR1 |
|--------------|--------------------------------------------------------------------------------------------------------------------------------|--|-----------|-------|-------|-------|-----|
| 281928       | <b>LoFlo sidestream CO2 sensor</b>                                                                                             |  | 1 pc      | X     | X     | X     |     |
| 282157       | <b>CAPNOSTAT®-5 CO2 sensor (90° angled)</b> . CO2 mainstream sensor for volumetric capnography                                 |  | 1 pc      | X     |       |       |     |
| 281718       | <b>CAPNOSTAT®-5 CO2 sensor</b> . CO2 Mainstream sensor for volumetric capnography                                              |  | 1 pc      | X     | X     | X     |     |
| 8000.0363    | <b>CO2 sidestream adapter</b> . Length: 2 m. ET > 4 mm. Dead space: 7 cc                                                       |  | Box of 10 | X     | X     | X     |     |
| 8000.0364    | <b>CO2 sidestream adapter. Infant/Neo</b> . ET ≤ 4 mm. Dead space: 1 cc                                                        |  | Box of 10 | X     | X     | X     |     |
| 8000.0260.01 | <b>CO2 Mainstream airway adapter</b> . ET > 4 mm. Dead space: 5 cc                                                             |  | Box of 10 | X     | X     | X     |     |
| 8000.0261.01 | <b>CO2 Mainstream airway adapter. Infant/Neo</b> . ET ≤ 4 mm. Dead space: < 1 cc. <i>Requires Adapter OD15 x ID15 (281803)</i> |  | Box of 10 | X     | X     | X     |     |
| 8000.0262.01 | <b>CO2 Mainstream airway adapter</b> . ET > 4 mm. Dead space: 5 cc                                                             |  | 1 pc      | X     | X     | X     |     |
| 281803       | <b>Adapter OD15/ID15</b> . Connect neonatal airway adapter to neonatal flow sensor                                             |  | Box of 25 | X     | X     | X     |     |



Single Use



Reusable



Autoclavable



P203A



P225G



P225F



P310C



950200



## SpO2

|        |                                                                             |  |           | C1/T1 | C3/C6 | G5/S1 | MR1 |
|--------|-----------------------------------------------------------------------------|--|-----------|-------|-------|-------|-----|
| P203A  | <b>SpO2 adult finger probe.</b> Single use. Length: 80 cm                   |  | Box of 24 | X     | X     | X     |     |
| P225G  | <b>SpO2 multi-site probe.</b> Reusable. Length: 160 cm                      |  | 1 pc      | X     | X     | X     |     |
| P225F  | <b>SpO2 finger probe.</b> Reusable. Length: 160 cm                          |  | 1 pc      | X     | X     | X     |     |
| P310C  | <b>SpO2 fingertip probe.</b> For patients > 50 kg. Reusable. Length: 160 cm |  | 1 pc      | X     | X     | X     |     |
| 950200 | <b>Cable holder</b>                                                         |  | 1 pc      | X     | X     | X     |     |



Single Use



Reusable



Autoclavable



1320



1330



1340



1350



1380 / 1381



10076606



10076605



10076604



10076608



## Nuflow - High Flow Neo


|      |                                              |  |           |   |   |   |   |
|------|----------------------------------------------|--|-----------|---|---|---|---|
| 1320 | <b>Nasal Cannula Set.</b> Small. Disposable  |  | Box of 10 | X | X | X | X |
| 1330 | <b>Nasal Cannula Set.</b> Medium. Disposable |  | Box of 10 | X | X | X | X |
| 1340 | <b>Nasal Cannula Set.</b> Large. Disposable  |  | Box of 10 | X | X | X | X |
| 1350 | <b>Nasal Cannula Set.</b> Xlarge. Disposable |  | Box of 10 | X | X | X | X |
| 1380 | <b>Fixation Patch.</b> Small. Disposable     |  | 1 pc      | X | X | X | X |
| 1381 | <b>Fixation Patch.</b> Large. Disposable     |  | 1 pc      | X | X | X | X |

## In2Flow - High Flow Pediatric/Adult

|          |                                     |  |           | C1/T1 | C2/C3/C6 | G5/S1 | MR1 |
|----------|-------------------------------------|--|-----------|-------|----------|-------|-----|
| 10076606 | <b>In2Flow Canulla.</b> Small       |  | Box of 20 | X     | X        | X     | X   |
| 10076605 | <b>In2Flow Canulla.</b> Medium      |  | Box of 20 | X     | X        | X     | X   |
| 10076604 | <b>In2Flow Canulla.</b> Large       |  | Box of 20 | X     | X        | X     | X   |
| 10076608 | <b>In2Flow Canulla</b> Tracheostomy |  | Box of 20 | X     | X        | X     | X   |







230185



230154



282016



161435 / 161436



260174 / 260176

## Transpulmonary pressure

|        |                                                                                                                                                                                  |   |   |          | C1/T1/MR1 | C6 | G5/S1 | IntelliCuff |
|--------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---|---|----------|-----------|----|-------|-------------|
| 230185 | <b>Esophageal catheter.</b> Set contains a catheter of 86 cm in length with a closed end and a balloon 9.5 cm in length                                                          | ⊗ | ♂ | Box of 5 |           | X  | X     |             |
| 230154 | <b>NutriVent multifunction naso-gastric catheter.</b> Enteral feeding tube with one balloon for esophageal pressure measurement with Hamilton Medical ventilator connecting line | ⊗ | ♂ | Box of 5 |           | X  | X     |             |

## Intellicuff

|        |                                                                                            |   |   |           |  |   |   |   |
|--------|--------------------------------------------------------------------------------------------|---|---|-----------|--|---|---|---|
| 282016 | <b>Cuff pressure tube for IntelliCuff 200 cm.</b> With bacterial filter and shut-off valve | ⊗ | ♂ | Box of 10 |  | X | X | X |
|--------|--------------------------------------------------------------------------------------------|---|---|-----------|--|---|---|---|

## Protective sleeve

|        |                                                                                        |   |      |   |  |  |  |  |
|--------|----------------------------------------------------------------------------------------|---|------|---|--|--|--|--|
| 161435 | <b>Ventilation hose protective sleeve 170 cm.</b> Washable cover for breathing circuit | ♂ | 1 pc | X |  |  |  |  |
| 161436 | <b>Ventilation hose protective sleeve 230 cm.</b> Washable cover for breathing circuit | ♂ | 1 pc | X |  |  |  |  |

## nCPAP pressure line

|        |                                                         |   |   |           |   |  |  |  |
|--------|---------------------------------------------------------|---|---|-----------|---|--|--|--|
| 260174 | <b>Pressure line for nCPAP/nCPAP-PC.</b> Length: 160 cm | ⊗ | ♂ | Box of 10 | X |  |  |  |
| 260176 | <b>Pressure line for nCPAP/nCPAP-PC.</b> Length: 310 cm | ⊗ | ♂ | Box of 10 | X |  |  |  |



Single Use



Reusable



Autoclavable