

## ZOLL accessories for R Series, X Series & Propaq MD



### CONTENT

- 2 | Electrodes
- 3 | Electrodes pediatric
- 3 | Electrodes with Real CPR Help®
- 4 | ECG electrodes and ECG accessories
- 5 | Internal handles with integrated paddles (spoons)
- 5 | External paddles
- 6 | Overmolded electrode paddles (spoons)
- 6 | Cables and connectors
- 7 | SurePower Battery System
- 8 | NIBP accessories
- 9 | Invasive Pressure accessories
- 9 | Temperature accessories
- 10 | ETCO2 accessories
- 11 | SpO2 accessories (1)
- 12 | SpO2 accessories (2)
- 12 | SpO2 accessories with Rainbow, SpO2, SpCO, SpMet
- 13 | Mounting solutions and storage
- 14 | X Serie Wall Mount
- 14 | X Serie Pole Mount for Airplane
- 15 | Bags for transport
- 15 | Other accessories
- 16 | Training accessories



8900.0211.01



8900.0212.01 / 8900.0222.01



8900.4004 / 8900.4003.49




8900.0801.01 / 8900.0802.01



8900.2105.01

## Electrodes

			R Series	X Series & Propaq MD
8900.0211.01	OneStep Basic Resuscitation Electrode. Multi-Function. Adult	box of 8	X	
8900.0212.01	OneStep Pacing Resuscitation Electrode. Multi-Function. Adult	box of 8	X	
8900.0222.01	OneStep Pacing Resuscitation Electrode. Multi-Function. Adult	1 pc	X	
8900.4004	Stat-Padz Electrode. Multi-Function. Adult	1 pc	X	X
8900.4003.49	Stat-Padz Electrode. Multi-Function. Adult	box of 12	X	X
8900.5002	Stat-Padz Electrode. Defib Only. Adult	box of 12	X	X
8900.0801.01	Stat-Padz II Electrode. Multi-Function. Adult	1 pc	X	X
8900.0802.01	Stat-Padz II Electrode. Multi-Function. Adult	box of 12	X	X
8900.2105.01	Pro-Padz Liquid Gel Radiolucent Electrode. Multi-Function. Adult	box of 12	X	X
8900.2303.01	Pro-Padz Biphasic Electrode. Multi-Function. Adult	box of 12	X	X
8900.2101.01	Pro-Padz Cardiology Speciality. Multi-Function. Adult	1 pc	X	X
8900.2100.01	Pro-Padz Cardiology Speciality. Multi-Function. Adult	box of 12	X	X
8900.4005	Pro-Padz Solid Gel Radiolucent. Multi-Function. Adult	box of 12	X	X





8900.0002.1901 / 8900.0002.2001



8900.0215.40 / 8900.0218.40



8900.0810.01



8900.3001.01 / 8900.3000.49



8900.0214.01



8900.0217.01 / 8900.0225.01



8900.0400 / 8900.0402



8900.0800.01

## Electrodes pediatric

			R Series	X Series & Propaq MD
8900.0002.1901	OneStep Pediatric CPR Electrode. Pediatric	1 pc	X	X
8900.0002.2001	OneStep Pediatric CPR Electrode. Pediatric	box of 8	X	X
8900.0215.40	OneStep Pediatric Resuscitation Electrode. Pediatric	box of 8	X	X
8900.0218.40	OneStep Pediatric Resuscitation Electrode. Pediatric	1 pc	X	X
8900.0810.01	Pedi-Padz II Electrode. Pediatric. 2 years shelf-life	1 pc	X	X
8900.3001.01	Pedi-Padz Electrode. Solid Gel. Multi-Function. Pediatric	1 pc	X	X
8900.3000.49	Pedi-Padz Electrode. Solid Gel. Multi-Function. Pediatric	box of 6	X	X

## Electrodes with Real CPR Help®

			R Series	X Series & Propaq MD
8900.0214.01	OneStep CPR complete electrode. Multi-Function. Adult	box of 8	X	
8900.0217.01	OneStep CPR A/A Electrode. Adult	box of 8	X	
8900.0225.01	OneStep CPR A/A Electrode. Adult	1 pc	X	
8900.0400	CPR Stat-Padz Electrode. Adult. Requires adapter 8000.0370	box of 8		X
8900.0402	CPR Stat-Padz Electrode. Adult. Requires adapter 8000.0370	1 pc		X
8900.0800.01	CPR-D Padz one piece defibrillation and CPR System. Requires adapter 8000.0370	1 pc	X	X



8900.0003 / 8900.0003.27



8900.0004



8900.1003.01




8900.0005



8900.0006

## ECG electrodes and ECG accessories

			R Series	X Series & Propaq MD
8000.0026	3-lead ECG Patient Cable with integral lead wire. 6 ft.	1 pc	X	
8009.0762.12	3-lead ECG. Cable Adapter. IEC	1 pc	X	
8000.1005.01	5-lead Patient cable with integral lead wires	1 pc	X	
8300.0800.12	Replacement 3-Lead cable. IEC	1 pc		X
8300.0801.12	5-Lead ECG Cable. IEC	1 pc		X
8300.0802.12	12-Lead one step ECG Cable. IEC includes 4-Lead trunk cable and removeable precordial 6 lead set	1 pc		X
8900.0003	3 ECG electrodes/pouch. 200 pouches	600 electrodes	X	X
8900.0004	4 ECG electrodes/pouch. 120 pouches	480 electrodes	X	X
8900.0005	5 ECG electrodes/pouch. 100 pouches	500 electrodes	X	X
8900.0006	6 ECG electrodes/pouch. 100 pouches	600 electrodes	X	X
8900.1003.01	Pediatric 3 ECG Electrode. 100 pouches	300 electrodes	X	X
8900.0704	Radiolucent 3 ECG Electrode (1,5" diameter). 30 pouches	300 electrodes	X	





Internal handles with integrated paddles (spoons)



8011.0500



8011.0501



External paddles

### Internal handles with integrated paddles (spoons)

			R Series	X Series & Propaq MD
8011.0140.01	Sterilizable Internal Handle. Without Switch. 1.0" spoon	1 pc	X	X
8011.0140.02	Sterilizable Internal Handle. Without Switch. 1.6" spoon	1 pc	X	X
8011.0140.03	Sterilizable Internal Handle. Without Switch. 2.7" spoon	1 pc	X	X
8011.0140.04	Sterilizable Internal Handle. Without Switch. 3.0" spoon	1 pc	X	X
8011.0140.05	Sterilizable Internal Handle. Without Switch. 2.0" spoon	1 pc	X	X
8011.0500	Autoclavable Internal Handles with No Button	1 pc	X	
8011.0501	Autoclavable Internal Handles with Button	1 pc	X	

### External paddles

8000.1010.XX	External Paddles w/controls and built-in pediatric electrodes	1 pc	X	X
--------------	---	------	---	---



Overmolded electrode paddles (spoons)



8000.0009.0301



8000.0370



80009.0020

## Overmolded electrode paddles (spoons)

			R Series	X Series & Propaq MD
8000.0277.01	2.7" (6.8 cm) dia. by 5.5" (14.3 cm). Short pair. Adult	1 pc	X	
8000.0277.02	2.7" (6.8 cm) dia. by 8.1" (20.6 cm). Long pair. Adult	1 pc	X	
8000.0277.03	2.0" (5.1 cm) dia. by 4.9" (12.5 cm). Short pair. Pediatric	1 pc	X	
8000.0277.04	2.0" (5.1 cm) dia. by 7.3" (18.7 cm). Long pair. Pediatric	1 pc	X	
8000.0277.05	1.6" (4.0 cm) dia. by 4.5" (11.4 cm). Short pair. Pediatric	1 pc	X	
8000.0277.06	1.6" (4.0 cm) dia. by 7.0" (17.7 cm). Long pair. Pediatric	1 pc	X	
8000.0277.07	1.0" (2.6 cm) dia. by 3.9" (10.0 cm). Short pair. Pediatric	1 pc	X	
8000.0277.08	1.0" (2.6 cm) dia. by 6.4" (16.3 cm). Long pair. Pediatric	1 pc	X	
8000.0277.09	3.0" (7.6 cm) dia. by 5.9" (15.0 cm). Short pair. Adult	1 pc	X	
8000.0277.10	3.0" (7.6 cm) dia. by 8.4" (21.3 cm). Long pair. Adult	1 pc	X	

## Cables and connectors

			R Series	X Series & Propaq MD
8000.0009.0301	X Series Pig Tail/AUX power break cable	1 pc		X
8000.0370	CPR Connector (Adapter)	1 pc		X
8009.0020	CPR-D Padz and CPR Stat Padz connector (for R Series and M Series)	1 pc	X	
8009.0749	OneStep CPR Cable (Support Real CPR Help)	1 pc	X	
8009.0750	OneStep pacing cable (Support both Real CPR Help and OneStep pacing)	1 pc	X	
8300.0006.76	X Series OneStep cable for OneStep electrodes	1 pc		X
8300.0783.01	X Series Multifunction Therapy cable. Defib/pacing cable	1 pc		X





8000.0580.01



8019.0535.01



8050.0030.XX



8300.0250.01



8200.0001.00XX

## SurePower Battery System

			R Series	X Series & Propaq MD
8000.0580.01	SurePower II rechargeable 6 hours battery	1 pc		X
8019.0535.01	SurePower rechargeable lithium ion battery pack	1 pc	X	
8050.0030.XX	SurePower 4 bay charging system (requires 4 x 8300.0250.01 to charge SurePower battery)	1 pc	X	X
8200.0001.00XX	SurePower single bay charger	1 pc	X	X
8300.0004	Replacement AC Power adaptor/charger	1 pc		X
8300.0250.01	SurePower charger adapter for SurePower II	1 pc		X
8700.0511.XX	SurePower Manager (Maintenance version)	1 pc	X	X



8000.0662



ZOLL Cuffs



REUSE Cuffs



SOFT Cuffs

## NIBP accessories

			R Series	X Series & Propaq MD
8000.0655	Hose. Air. 1.5 m. NIBP	1 pc	X	
8000.0662	Hose. Air. 3 m. NIBP	1 pc	X	
8000.1650	Cuff. ZOLL. All purpose. Pediatric/Small Adult. 17 - 25 cm	1 pc	X	
8000.1651	Cuff. ZOLL. All purpose. Adult. 23 - 33 cm	1 pc	X	
8000.1653	Cuff. ZOLL. All purpose. Large Adult. 31 - 40 cm	1 pc	X	
8000.1654	Cuff. ZOLL. All purpose. Thigh. 38 - 50 cm	1 pc	X	
8000.1655	Cuff. ZOLL. All purpose. Child. 12 - 19 cm	1 pc	X	
8300.0002.01	Hose. Dual lumen. NIBP. 3 m	1 pc		X
8300.0002.02	Hose. Dual lumen. NIBP. 1,5 m	1 pc		X
8300.0796.01	Hose. Infant/neonate. 2,4 m. Single lumen NIBP. W/female luer cuff connector	1 pc		X
8300.0797.01	Neonatal Cuff kit. One each sizes #1 - #5. Single tube w/male luer connector	1 pc		X
REUSE-07-2MQ (WA)	Cuff. REUSE. Infant. 2-Tube. Twist lock connector. 9-13 cm	1 pc		X
REUSE-08-2MQ (WA)	Cuff. REUSE. Small Child. 2-Tube. Twist lock connector. 12-16 cm	1 pc		X
REUSE-09-2MQ (WA)	Cuff. REUSE. Child. 2-Tube. Twist lock connector. 15-21 cm	1 pc		X
REUSE-10-2MQ (WA)	Cuff. REUSE. Small Adult. 2-Tube. Twist lock connector. 20-26 cm	1 pc		X
REUSE-11-2MQ (WA)	Cuff. REUSE. Adult. 2-Tube. Twist lock connector. 25-34 cm	1 pc		X
REUSE-11L-2MQ (WA)	Cuff. REUSE. Adult Long. 2-Tube. Twist lock connector. 25-34 cm	1 pc		X
REUSE-12-2MQ (WA)	Cuff. REUSE. Large Adult. 2-Tube. Twist lock connector. 32-43 cm	1 pc		X
REUSE-12L-2MQ (WA)	Cuff. REUSE. Large Adult Long. 2-Tube. Twist lock connector. 32-43 cm	1 pc		X
REUSE-13-2MQ (WA)	Cuff. REUSE. Thigh. 2-Tube. Twist lock connector. 40-55 cm	1 pc		X
SOFT-07-2MQ (WA)	Disposable Cuff. SOFT. Infant. 2-Tube. Twist lock connector. 9-13 cm	box of 20		X
SOFT-08-2MQ (WA)	Disposable Cuff. SOFT. Small child. 2-Tube. Twist lock connector. 12-16 cm	box of 20		X
SOFT-09-2MQ (WA)	Disposable Cuff. SOFT. Child. 2-Tube. Twist lock connector. 15-21 cm	box of 20		X
SOFT-10-2MQ (WA)	Disposable Cuff. SOFT. Small adult. 2-Tube. Twist lock connector. 20-26 cm	box of 20		X
SOFT-11-2MQ (WA)	Disposable Cuff. SOFT. Adult. 2-Tube. Twist lock connector. 25-34 cm	box of 20		X
SOFT-11L-2MQ (WA)	Disposable Cuff. SOFT. Adult Long. 2-Tube. Twist lock connector. 25-34 cm	box of 20		X
SOFT-12-2MQ (WA)	Disposable Cuff. SOFT. Large Adult. 2-Tube. Twist lock connector. 32-43 cm	box of 20		X
SOFT-12L-2MQ (WA)	Disposable Cuff. SOFT. Large Adult long. 2-Tube. Twist lock connector. 32-43 cm	box of 20		X
SOFT-13-2MQ (WA)	Disposable Cuff. SOFT. Thigh. 2-Tube. Twist lock connector. 40-55 cm	box of 20		X







6840.85



8000.0668



8000.0669



8000.0670



8000.0672



8000.0674

## Invasive Pressure accessories

			R Series	X Series & Propaq MD
8300.0787.01	Transducer interface cable with right angle connector. Edwards	1 pc		X
8300.0788.01	Transducer interface cable with right angle connector. Abbott	1 pc		X
6840.85	Invasive blood pressure cable BD. Argon Medical	1 pc		X

## Temperature accessories

8000.0668	YSI reusable. Adult. Esophageal Rectal Temperature probe	1 pc		X
8000.0669	YSI reusable. Pediatric. Esophageal Rectal probe	1 pc		X
8000.0670	YSI reusable. Adult. Skin Temperature probe	1 pc		X
8000.0671	YSI reusable. Pediatric. Skin Temperature probe	1 pc		X
8000.0672	YSI single use. Adult. Esophageal Rectal probe	1 pc		X
8000.0673	YSI single use. Adult. Skin Temperature probe	1 pc		X
8000.0674	Disposable Temperature Sensor Adapter Cable	1 pc		X




8000.0260.01



8000.0262.01



8000.0312



## ETCO2 accessories

			R Series	X Series & Propaq MD
XS04620	Filterline Set. Adult/Pediatric	box of 25		X
XS04624	Filterline H Set. Adult/Pediatric	box of 25		X
006324	Filterline H Set. Infant/Neonatal	box of 25		X
010787	Vitaline H Set. Adult/Pediatric	box of 25		X
009822	Smart CapnoLine Plus O2 Adult (O2 tubing)	box of 25		X
007269	Smart CapnoLine Plus O2 Pediatric (O2 tubing)	box of 25		X
009818	Smart CapnoLine Plus (O2 connector)	box of 25		X
007266	Smart CapnoLine Pediatric	box of 25		X
8000.0260.01	Single Patient Use Pediatric/Adult Airway Adapter	box of 10	X	
8000.0261.01	Single Patient Use Neonatal/Pediatric Airway Adapter	box of 10	X	
8000.0262.01	Reusable Adult Airway Adapter	1 pc	X	
8000.0263.01	Reusable Neonatal Airway Adapter	1 pc	X	
8000.0265.01	Single Patient Use Pediatric Adult Airway Adapter with Mouthpiece	box of 10	X	
8000.0266.01	Sensor Cable Holding Clips	box of 5	X	
8000.0312	CAPNO 5 Mainstream CO2 Sensor and Cable	1 pc	X	



8000.0004.56



8000.0293 / 8000.0298



8000.0294



8000.0295



8000.0320



8000.0321



8000.0322



8000.0323



8000.0324



8000.0325

## SpO2 accessories (1)

			R Series	X Series & Propaq MD
16-1895	EarClip SpO2 Sensor. Reusable. RED	1 pc	X	X
8000.0004.56	E1 Ear Single Patient Sensor. RED LNCS Sensor. Masimo	Box of 10	X	X
8000.0004.59	M-LNCS DCI Reusable Sensor	1 pc		X
8000.0293	SpO2 LNCS Reusable Patient Cable. 3 m	1 pc	X	
8000.0294	SpO2 LNCS Adult Reusable Sensor	1 pc	X	X
8000.0295	SpO2 LNCS Pediatric Reusable Sensor (10 < 50 kg)	1 pc	X	X
8000.0296	SpO2 LNCS Adult Reusable One-Piece Sensor and Cable. 3,6 m	1 pc	X	
8000.0298	SpO2 LNCS Reusable Patient Cable. 1,2 m	1 pc	X	
8000.0320	SpO2 LNCS Disposable Adult Sensors (> 30 kg)	box of 20	X	
8000.0321	SpO2 LNCS Disposable Pediatric Sensors (10 < 50 kg)	box of 20	X	
8000.0322	SpO2 LNCS Disposable Infant Sensors (3 < 20 kg)	box of 20	X	
8000.0323	SpO2 LNCS Disposable Neonatal Sensors (< 3 kg)	box of 20	X	
8000.0324	SpO2 LNCS Disposable Preterm Neonatal Sensors	box of 20	X	
8000.0325	SpO2 LNCS Extension Cable with DB-9 Connector	1 pc	X	



8000.0326



8000.0327



8000.0004.57



8000.0330 / 8000.0331



8000.0332 / 8000.0334



8000.0333 / 8000.0335



8000.0343 / 8000.0344



8000.0345 / 8000.0346

## SpO2 accessories (2)

			R Series	X Series & Propaq MD
8000.0326	SpO2 LNOP Sensor to LNCS Adapter Cable	1 pc		X
8000.0327	SpO2 LNCS Sensor to LNOP Adapter Cable	1 pc		X
8000.0330	Red. 1,2 m. Reusable Patient Cable. LNC-04	1 pc		X
8000.0331	Red. 3,0 m. Reusable Patient Cable. LNC-10	1 pc		X
8000.0332	Red. 0,9 m. DCI. Adult. Reusable Patient Cable/Sensor	1 pc		X
8000.0333	Red. 0,9 m. DCI. Pediatric. Reusable Patient Cable/Sensor	1 pc		X
8000.0334	Red. 3,6 m. DCI. Adult. Reusable Patient Cable/Sensor	1 pc		X
8000.0335	Red. 3,6 m. DCI. Pediatric. Reusable Patient Cable/Sensor	1 pc		X

## SpO2 accessories with Rainbow, SpO2, SpCO, SpMet

8000.0003.71	Rainbow. Adult. Reusable Patient Sensor for SpO2, SpCO, SpMet	1 pc		X
8000.0003.72	Rainbow. Pediatric. Reusable Patient Sensor for SpO2, SpCO, SpMet	1 pc		X
8000.0004.57	E1 Ear Single Patient Sensor. Rainbow M-LNCS. Masimo	Box of 10	X	X
8000.0341	Rainbow RC-4. 4' (1,2m). Reusable Patient Cable	1 pc		X
8000.0342	Rainbow RC-12. 12' (3,6 m). Reusable Patient Cable	1 pc		X
8000.0343	Rainbow. 8' (2,4 m). DCI. Adult. Reusable Patient Cable/Sensor	1 pc		X
8000.0344	Rainbow. 12' (3,6 m). DCI. Adult. Reusable Patient Cable/Sensor	1 pc		X
8000.0345	Rainbow. 8' (2,4 m). DCI. Pediatric. Reusable Patient Cable/Sensor	1 pc		X
8000.0346	Rainbow. 12' (3,6 m). DCI. Pediatric. Reusable Patient Cable/Sensor	1 pc		X



0006466



1041.1000.10



1041.1000.11



1041.1000.20



21041.1000.21



8000.0979.01



8000.0980.01




8000.0981.01

## Mounting solutions and storage

			R Series	X Series & Propaq MD
0004066	Swivel plate for 0006466 X-Series Bracket without power connection. Profile Vehicles	1 pc		X
0006466	X-Series Bracket without power connection. EN 1789 certified. Profile Vehicles	1 pc		X
1041.1000.10	X Series Bracket Wall Mount without power connection. EN 1789 certified. CRS	1 pc		X
1041.1000.11	X Series Bracket Wall Mount with power connection. EN 1789 certified. CRS	1 pc		X
1041.1000.20	X Series Bracket Desktop Mount with power connector. EN 1789 certified. CRS	1 pc		X
1041.1000.21	X Series Bracket Desktop Mount without power connector. EN 1789 certified. CRS	1 pc		X
8000.0975.01	R Series Channel Wall Mount 19"	1 pc	X	
8000.0976.01	R Series Channel Wall Mount 31"	1 pc	X	
8000.0977.01	R Series Roll Stand. Includes Rolling Mount, Handle, and Accessory Basket	1 pc	X	
8000.0979.01	R Series Counter Top Mount. Used to Connect the R Series to a counter top	1 pc	X	
8000.0980.01	R Series Wall Arm Mount. Used to connect the R Series to the Channel Wall Mount	1 pc	X	
8000.0981.01	R Series Base Plate. Allows for use with the GCX Mounting System	1 pc	X	




FM300A100



FM100E100-7



FM300A102-3



FM300A103

### X Serie Wall Mount (All 3 parts below must be ordered at same time)

			R Series	X Series & Propaq MD
FM100E100-7	Mini Ferno Track 7"	1 pc		X
FM300A100	Mini ProV Track Mount (Mini ProV Track Receiver)	1 pc		X
FM300A102-3	Mini ProV Adapter (Propaq MD/X-serie)	1 pc		X

### X Serie Pole Mount for Airplane (All 3 parts below must be ordered at same time)

FM100A103	Mini Ferno Track Pole Mount (Ask for diameter before ordering)	1 pc		X
FM300A100	Mini ProV Track Mount (Mini ProV Track Receiver)	1 pc		X
FM300A102-3	Mini ProV Adapter (Propaq MD/X-serie)	1 pc		X





1042.1010.05



1042.1010.04



8000.0004.6301



8000.0963



8000.0964



8707.0005.0201



8000.0020.0501

## Bags for transport

			R Series	X Series & Propaq MD
1042.1010.05	X Series Softbag. Expandable Black - Shoulder Strap must be ordered seperatly (item# 1042.1010.04)	1 pc		X
1042.1010.04	X Series Shoulder Strap set	1 pc		X
8000.0004.6301	R Series Transport Pack/bag	1 pc	X	
8000.0963	R Series Electrode Pouch (for right sided attachment)	1 pc	X	
8000.0964	R Series Transport Pouch (for left sided attachment)	1 pc	X	
8000.0984.01	Left Side Electrode Pocket for R Series Transport Pack	1 pc	X	
8000.0985.01	Right Side Accessories Pocket for R Series Transport Pack	1 pc	X	
8000.0986.01	Top Accessories for R Series Transport Pack	1 pc	X	
8000.0987.01	Bed Rail Hooks with Shoulder Strap for R Series Transport Pack	1 pc	X	
8707.0005.0201	X Series Carry Case. ZOLL	1 pc		X

## Other accessories

8000.0008.7501	ECG plain white paper. 80 mm	box of 6		X
8000.0008.7701	Recorder Paper. 90 mm. Fan Fold	10 packages	X	
8000.0020.0501	Cable Sleeve. Royal Blue	1 pc		X
9143.0001	128 MB Compact Flash Data Card (for use with R Series only)	1 pc	X	
9214.0203	Ambicom WL1100C-CF 802.11b WiFi Card	1 pc	X	



8009.0751.01



8012.0206



8900.0008.6101



8900.5005



8900.0190 / 8900.0195



8900.0805.01



8900.5007

## Training accessories

			R Series	X Series & Propaq MID
8009.0751.01	SeeThru CPR Simulator	1 pc	X	X
8012.0206	12-Lead Simulator	1 pc		X
8900.0008.6101	Training Pedi-Padz II Electrode. Pediatric. 2 years shelf-life. Requires adapter	box of 6	X	X
8900.0008.6201	One Step Pediatric CPR Training Kit	1 pc	X	
8900.0190	Training CPR Stat-Padz Electrode. Y connector for simulator connection. Requires adapter	1 pc	X	X
8900.0195	Replacement CPR Stat-Padz Electrodes	box of 8	X	X
8900.0240.01	OneStep CPR Training Cable and Electrode (CPR A/A)	1 pc	X	
8900.0245.01	Replacement OneStep Training Electrode (CPR A/A) for 8900.0240.01	box of 8	X	
8900.0805.01	Training Stat-Padz II Electrode. Adult. 2 years shelf-life. Requires adapter	box of 6	X	X
8900.5005	Stat-Padz training electrodes. Multi-Function	box of 6		X
8900.5007	CPR-D-Padz Demo Pad. Adult. Green AED Connector	1 pc		X

