


X SERIES® ADVANCED

The new best-in-class monitor/
defibrillator from ZOLL®



ZOLL's new X Series® Advanced monitor/defibrillator helps providers manage patients more effectively than ever before. X Series Advanced is the first all-in-one device of its kind to provide the clinical data EMS teams need at every step to improve outcomes.



Real BVM Help™ and TBI Dashboard

With groundbreaking new technology available only from ZOLL, the X Series Advanced finally addresses some of the key challenges providers face in both pre-hospital and hospital environments today.



Real BVM Help – Real-time feedback for ventilation

With Real BVM Help on the X Series Advanced, providers can view delivered tidal volume and rate directly on the monitor. This real-time feedback guides users to deliver high-quality ventilation in a manner that has never been possible until now.

Real BVM Help on the ZOLL X Series Advanced gives providers:

- Target and delivered ventilation volume
- Target and delivered ventilation rate
- A ventilation quality indicator and countdown timer



TBI Dashboard: The information you need when it matters most


Traumatic Brain Injury (TBI) patients present a unique set of challenges. Early treatment is critical, as additional secondary injury often occurs after the initial event.¹

ZOLL's TBI Dashboard on the X Series Advanced monitor/defibrillator provides trending information on the parameters most critical to a TBI patient, enabling clinicians to quickly detect potential patient deterioration.

¹Spaite DW, et al. Acad Emerg Med. 2014;21:7:818-830


Meet your telemedicine needs with **Remote View** on the X Series Advanced

With Remote View, clinicians can view data from the X Series Advanced from a remote location. This telemedicine solution from ZOLL enables effective decision support to allow EMS providers meet their telehealth needs.



Real BVM Help integrates with existing ZOLL CPR technologies

For cardiac arrest patients, the X Series Advanced monitor/defibrillator offers both Real CPR Help® and Real BVM Help, giving providers simultaneous real-time feedback on both CPR and ventilation quality. When used with adult ZOLL CPR electrodes, the X Series Advanced device displays the Real BVM Help Dashboard next to the Real CPR Help Dashboard.



Improve your system's performance with **RescueNet CaseReview**

RescueNet® CaseReview gives providers the ability to view data from the X Series Advanced post-case for effective QA/QI. Data available in CaseReview includes ventilation performance with Real BVM Help, CPR performance with Real CPR Help®, and automatic recognition and integration of performance with ZOLL's automated CPR device, the AutoPulse®.

X Series Advanced: The new best-in-class monitor/defibrillator from ZOLL.

For more information, visit zoll.com/xseriesadvanced